



Rosi


Quick Guide

WATCH THE
TUTORIAL VIDEO
HERE.



01. Charging

Charge the device fully before first use. Connect the USB Type-C cable to the device, then plug the other end into a standard AC adapter or a computer USB port. Note: An AC adapter is not included. Please use a commercially available adapter.




02. Setup

Connect the butterfly pad (electrode pad) to the device. A clicking sound will indicate it is properly connected.



03. Attachment

Peel off the protective film on the back of the butterfly pad and attach it to the area of skin where you feel pain. Do not discard the protective film, as it will be used for storage.




04. Rosi Power ON

Press the center ON/OFF button to turn on the device. The indicator light will turn green and a beep sound will be heard.









05. Mode Selection & Intensity Adjustment


After turning on the device, press the ON/OFF button to switch modes. Mode 1 operates for 30 minutes, Modes 2 and 3 operate for 60 minutes continuously. Mode 1: 1 flash, Mode 2: 2 flashes, Mode 3: 3 flashes. You can adjust the intensity using the "+" and "-" buttons.



Mode Selection and Adjustment

Setting Modes ()					
Step	Mode	Operation	Indicator	Sound (Short Beep)	Operation Time
①	On	Press 	Steady green light	Once	30min
②	Mode 1				
③	Mode 2	Press  again	Blinks twice, then steady green	Once	60min
④	Mode 3	Press  again	Blinks three times, then steady green	Once	60min
⑤	Mode 1	Press  again	Blinks once, then steady green	Once	30min
<ul style="list-style-type: none">• The device starts in Mode 1 upon powering on.• Pressing the  button cycles through Modes 1 to 3 sequentially.• Each mode stops automatically after the preset time ends.					

Output Level Adjustment(+/-)			
Action	Button	Sound	Range
To increase intensity	+	Short beep	Up to 30
To decrease intensity	-	No sound	Down to 30
Important Notes: <ul style="list-style-type: none">• Select the mode before adjusting intensity.• Each time the power is turned on, the output level resets to "0".• If powered on without being worn, the device will shut off automatically after 60 seconds.• Turn on the power after wearing the device.• If the + button is not pressed within 60 seconds of powering on, the device will shut off automatically.			

Turn Off			
Turn Off	 Press 2 sec	Green light off	No sound

Precautions for Use

[Important]

Please be sure to read "1. Safety Precautions (Contraindications, Prohibited Use, and Warnings)" in the instruction manual before use.

- Keep out of reach of children.
- Do not use if you are pregnant or may be pregnant.
- If you experience any symptoms such as rash, redness, or itching while using this product, stop using it immediately.
- Do not use in water or humid environments (such as saunas or hot bath facilities).
- Avoid using the product in high-altitude areas over 3,000 meters.
- Do not use while driving or operating machinery, or during any activity where unexpected muscle contractions could pose a risk.
- Do not share the butterfly pad with others.
- Do not use while sleeping.
- Ensure that no foreign objects enter the device.
- Do not use the device while it is charging.

Instruction Manual

For detailed instructions and important safety information, including "1. Safety Precautions (Contraindications, Prohibited Use, and Warnings)," please refer to the instruction manual. You can download it from the link below.



BIC GROUP Ltd
Shin-Kokusai Bldg 4F, 3-4-1 Marunouchi, Chiyoda-ku,
Tokyo 100-0005 Japan
Rosicare Customer Support
Email: contact@rosicare.com
☎ 03-6269-9120